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YIN Ge, FAN Jun, HUANG Taohong. Analysis of 18 polycyclic aromatic hydrocarbons in carbon black by GCMS isotope internal standard method [J]. Environmental Chemistry, 2021, 40(4): 1288-1291.

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Analysis of 18 polycyclic aromatic hydrocarbons in carbon black by GCMS isotope internal standard method

YIN Ge* FAN Jun HUANG Taohong

(Shimadzu China Co. Ltd, Shanghai, 200233, China)

Abstract: Polycyclic aromatic hydrocarbons (PAHs) are a class of organic pollutants with natural and anthropogenic sources. In this study, ultrasonic extraction coupled with isotope internal standard method was used to establish a fast and accurate GCMS method for the analysis of 18 PAHs in carbon black. InertCap-35 column was applied to obtain 18 PAHs components with good separation. The linearity of the standard curve is good in the concentration range of 5—1000 ng·mL⁻¹, and the correlation coefficients are all above 0.99. Standard solution at 10 ng·mL⁻¹ were injected continuously for 6 injections, and the peak area RSD was less than 8.79%, indicating good repeatability. The main components detected in carbon black samples were fluorene, fluoranthene and naphthalene. The method is simple, convenient, and can effectively analyze the contents of 18 polycyclic aromatic hydrocarbons in carbon black samples.

Keywords: PAH, GCMS-QP2020 NX, ultrasonic extraction, carbon black.

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Table 2 Correlation coefficients and limit of detection of PAHs

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Table 4 Results on PAHs in carbon black samples

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参考文献 (References)

- [1] LIU Y, BECKINGHAM A, RUEGNER T, et al. Comparison of sedimentary PAHs in the rivers of Ammer (Germany) and Liangtan (China) : Differences between early- and newly- industrialized countries [J]. Environ Sci Technol, 2013, 47 : 701-709
- [2] 田菲菲, 袁庆丹, 汪勇 等. 气相色谱-质谱法检测运动场地塑胶面层中的 18 种多环芳烃 [J]. 环境化学, 2016, 35 (12) : 2622-2625
- TIAN F F, YUAN Q D, WANG Y, et al. Determination of 18 PAHs in plastic surface of sport ground by GCMS [J]. Environmental Chemistry, 2016, 35 (12) : 2622-2625 (in Chinese).

bjhx.rcees.ac.cn